

Beijing Solarbio Science & Technology Co.,Ltd. One-stop solution for life science research.

# YPDA Medium (ready-to-use)

Cat No.: Y8063 Package: 500mL Storage: 2-8°C. Valid for at least one year

English name: Yeast Extract Peptone Dextrose Adenine Medium

## **Product brief introduction**

Composite medium sterile solution for yeast routine growth. It is mainly composed of yeast extract, peptone, glucose, etc., and does not contain AGAR. It is filtered and sterilized. Can be used for a variety of bacterial culture, mainly for yeast growth culture. YPDA medium is mainly composed of tryptone, yeast extract, glucose, adenine sulfate and so on.

### Product formula and usage:

Product ingredient	Concentration
Yeast Extract Powder	1%
Peptone	2%
Dextrose (glucose)	2%
Adenine sulfate	0.003%

Usage: Disinfect the bottle cap with 75% alcohol, and use it directly after opening the plastic cap under sterile conditions

#### This reagent is only used in scientific research.

### **Related products**

LA0223 YPD Liquid Medium (Ready-to-use)
Y8063 YPDA Liquid Medium (Ready-to-use)
LA0661 MD Liquid Medium Minimal Dextrose Medium
LA5252 BMGY/BMMY Medium (ready-to-use)

Note: For more information about this product, please refer to Solarbio website.

