

## Graphene oxide dispersion water, <500nm

**Cat:** NM000132


**Specification:** 20ml/50ml/100ml/500ml

**Storage:** Sealed, avoid light, and keep at 4°C. Expiry date: 6 months.

### Product Introduction:

Product name	Graphene oxide dispersion
Thickness	0.8-1.2nm
Diameter	< 500nm
Purity	~99%
Single layer ratio	~99%
Concentration	2mg/ml
Solvent	Water

### Representation Test Diagram:



TEM of graphene oxide water dispersion (<500 nm)

### Application:

Adsorption, Separation, Catalysis, Solid Phase Extraction, Detection and Sensors etc.

### Note:

1. This product is for scientific research use only. Do not use for medical, clinical diagnostic or therapeutic, food or cosmetic purposes. Do not store in ordinary residential areas.
2. For your safety and health, please wear a laboratory coat and wear disposable gloves and masks to operate.