

DMSA Coated Fe₃O₄ Nanoparticles

Cat: NM000610

Specification: 10mg/20mg


Storage: Sealed, avoid light, and keep at 4°C. Expiry date: 6 months.

Product Introduction:

Product name	DMSA Coated Fe ₃ O ₄ Nanoparticles
Size	~7, 20, 30, 50nm (TEM)
Color	Brown
Concentration	~4 mg/ml (Fe)

Note: The particle size is a single measurement result, and a small amount of deviation is allowed.

Representation Test Diagram:



TEM image of DMSA Coated Fe₃O₄ Nanoparticles (7nm)

Application:

Magnetic resonance imaging contrast, biomolecule coupling, immobilization, nano probe construction, etc.

Note:

1. This product is for scientific research use only. Do not use for medical, clinical diagnostic or therapeutic, food or cosmetic purposes. Do not store in ordinary residential areas.
2. For your safety and health, please wear a laboratory coat and wear disposable gloves and masks to operate.