

Ag Nanoparticles

Cat: NM000629


Specification: 10ml

Storage: Sealed, avoid light, and keep at 4°C. Expiry date: 6 months.

Product Introduction:

Product name	Ag Nanoparticles
Size	100nm
Concentration	0.1mg/ml
Ingredient	Ag Nanoparticles、 water

Representation Test Diagram:



UV absorption spectrum of Ag Nanoparticles (100nm)

Product Features:

1. Ag nanoparticle surface with a negative charge, which can be adjusted as needed
2. High specific surface area, high load, easy to functionalization on the surface
3. Good dispersibility and stability

Note:

1. This product is for scientific research use only. Do not use for medical, clinical diagnostic or therapeutic, food or cosmetic purposes. Do not store in ordinary residential areas.
2. For your safety and health, please wear a laboratory coat and wear disposable gloves and masks to operate.