

Au Nanoparticles 100nm

Cat: NM000669

Specification: 10ml

Storage: Sealed, avoid light, and keep at 4°C. Expiry date: 6 months.


Product Introduction:

Product name	Au Nanoparticles 100nm
Particle size	$\sim 100 \pm 20\text{nm}$
Concentration	0.05mg/ml
Maximum absorption wavelength	$\sim 560 \pm 10\text{nm}$


Representation Test Diagram:




SEM image of Au Nanoparticles (100nm)



TEM images of Au Nanoparticles (100nm)



TEM images of Au Nanoparticles (100nm)



UV - Visible absorption spectrum of Au Nanoparticles (100nm)

(Note: This spectrum is a single measurement data, allowing floating between different batches.)

Application:

Cell detection, microscopic tracer, immunoblotting, biosensor, biochip, drug residue detection, gene regulatory drug synthesis, drug transport, photochemical therapy etc.

Note:

1. This product is for scientific research use only. Do not use for medical, clinical diagnostic or therapeutic, food or cosmetic purposes. Do not store in ordinary residential areas.
2. For your safety and health, please wear a laboratory coat and wear disposable gloves and masks to operate.