

Nano gold colloid

Cat: NM000680


Specification: 50ml

Storage: Sealed, avoid light, and keep at 4°C. Expiry date: 6 months.


Product Introduction:

Product name	Nano gold colloid
Concentration	0.1g/L
Electron microscope size	10-20nm
Maximum absorption wavelength	518-521nm

Representation Test Diagram:



TEM of Nano gold colloid (10-20nm)



UV-Vis spectrum of Nano gold colloid (10-20nm)

Application:

Cell detection, microscopic tracer, immunoblotting, biosensor, biochip, drug residue detection, gene regulatory drug synthesis, drug transport, photochemical therapy etc.

Note:

1. This product is for scientific research use only. Do not use for medical, clinical diagnostic or therapeutic, food or cosmetic purposes. Do not store in ordinary residential areas.
2. For your safety and health, please wear a laboratory coat and wear disposable gloves and masks to operate.