

Nano Silicon Dioxide, KH570 Modified

Cat: NM000830


Specification: 100g/500g

Storage: Sealed, avoid light, and keep at room temperature. Expiry date: 1 year.


Product Introduction:

Product name	Nano silicon dioxide
Purity	99%
Average grain size	20nm
BET	140.21m ² /g
Appearance	White fluffy powder
Surface modification	KH570 modified

Representation Test Diagram:



SEM of Nano silicon dioxide (100nm)



TEM of Nano silicon dioxide (100nm)

Application:

UV resistant materials, textiles, photocatalytic catalysts, optical / magnetic memory, self-cleaning glass, sunscreens, paints, inks, food packaging materials, paper industry, aerospace industry, lithium batteries, solar cells and gas adsorbents etc.

Note:

1. This product is for scientific research use only. Do not use for medical, clinical diagnostic or therapeutic, food or cosmetic purposes. Do not store in ordinary residential areas.
2. For your safety and health, please wear a laboratory coat and wear disposable gloves and masks to operate.