

YPD Liquid Medium (Ready-to-use)

Cat No.: LA0223

Package: 500mL

Storage: Store at room temperature for at least one year.

English name: Yeast Extract Peptone Dextrose Medium

Product brief introduction

Composite medium sterile solution for yeast routine growth. It is mainly composed of yeast extract, peptone, glucose, etc., and does not contain AGAR. It is filtered and sterilized. YPD medium, also known as YEPD medium, can be used for a variety of bacterial cultures, mainly for the conventional growth of yeast composite medium. YPD medium is mainly composed of tryptone, yeast extract and glucose, and its concentration is higher than that of conventional LB medium, without AGAR.

Product formula

Product formula	Concentration
Yeast Extract Powder	1%
Peptone	2%
Dextrose (glucose)	2%

Usage: Disinfect the bottle cap with 75% alcohol, and use it directly after opening the plastic cap under sterile conditions

This reagent is only used in scientific research.

Related products

LA0223 YPD Liquid Medium (Ready-to-use)

Y8063 YPDA Liquid Medium (Ready-to-use)

LA0661 MD Liquid Medium Minimal Dextrose Medium

LA5252 BMGY/BMMY Medium (ready-to-use)

Note: For more information about this product, please refer to Solarbio website.

