

BMGY/BMMY Medium (ready-to-use)

Cat No.: LA5252

Package: 500mL

Storage: Store at 4-8°C for one year.

English name: BMGY/BMMY Medium

Product brief introduction

It is used to induce the expression of Pichia Pastoris, containing Biotin and non-amino yeast nitrogen (YNB).

Product formula and usage:

Product ingredient	Concentration
Yeast Extract Powder	1%
Peptone	2%
Amino free yeast nitrogen Source (YNB)	1.34%
Biotin (D-biotin)	0.4mg/L
Phosphate buffer system	100mM

It can be used to induce the expression of Pichia pastoris by adding 1% glycerol or 0.5% methanol.

Usage: Disinfect the bottle cap with 75% alcohol, and use it directly after opening the plastic cap under sterile conditions.

This reagent is only used in scientific research.

Related products

LA0223 YPD Liquid Medium (Ready-to-use)

Y8063 YPDA Liquid Medium (Ready-to-use)

LA0661 MD Liquid Medium Minimal Dextrose Medium

LA5252 BMGY/BMMY Medium (ready-to-use)

Note: For more information about this product, please refer to Solarbio website.

