

YPDA Medium (ready-to-use)

Cat No.: Y8063 Package: 500mL

Storage: 2-8°C. Valid for at least one year

English name: Yeast Extract Peptone Dextrose Adenine Medium

Product brief introduction

Composite medium sterile solution for yeast routine growth. It is mainly composed of yeast extract, peptone, glucose, etc., and does not contain AGAR. It is filtered and sterilized. Can be used for a variety of bacterial culture, mainly for yeast growth culture. YPDA medium is mainly composed of tryptone, yeast extract, glucose, adenine sulfate and so on.

Product formula and usage:

1100000101111010000000	
Product ingredient	Concentration
Yeast Extract Powder	1%
Peptone	2%
Dextrose (glucose)	2%
Adenine sulfate	0.003%

Usage: Disinfect the bottle cap with 75% alcohol, and use it directly after opening the plastic cap under sterile conditions

This reagent is only used in scientific research.

Related products

LA0223 YPD Liquid Medium (Ready-to-use) Y8063 YPDA Liquid Medium (Ready-to-use)

LA0661 MD Liquid Medium Minimal Dextrose Medium

LA5252 BMGY/BMMY Medium (ready-to-use)

Note: For more information about this product, please refer to Solarbio website.





Address: No. 85, Liandong U Valley, Middle Zone, 101102, Tongzhou Dist, Beijing, China