

## 卢戈碘液(标准碘液, 5%)

货号: G1069

规格: 100mL/500mL

保存: 2-8°C, 避光保存, 有效期 1 年。

### 产品介绍:

卢戈碘液又叫鲁哥氏染色液, 含 5% 的碘, 是碘和碘化钾的水溶液。常应用于检测叶片、食物等样品中的淀粉, 可将淀粉染成黑色或蓝黑色。用于细胞染色, 使细胞核着色从而更易观察。

### 相关产品:

G1066 卢戈碘液(0.33%, 革兰氏)

G1062 结晶紫染色液 (1%)

G1063 结晶紫染色液 (0.1%)

G1060 革兰氏染色试剂盒





## Lugol's Iodine Solution(5%, Standard Solution)

**Cat:** G1069

**Size:** 100mL/500mL

**Storage:** 2-8°C, avoid light, valid for 1 year.

### Introduction

Lugol's Iodine Solution(5%, Standard Solution) also known as Lugol's dye (5%), contains 5% iodine, is an aqueous solution of iodine and potassium iodide. It is often used in the detection of starch in leaves, food and other samples. The starch can be dyed black or blue-black. If used for cell staining, please stain the nucleus to make it easier to view.

### Related Products

G1066 *Lugol Solution, 0.33%, For Gram Stain*

G1062 *Crystal Violet Ammonium Oxalate Solution, 1%*

G1063 *Crystal Violet Ammonium Oxalate Solution, 0.1%*

G1060 *Gram Stain Kit*

