REFERENCES '

- 1. Page-McCaw, A. et al. (2007) Nat. Rev. Mol. Cell. Biol. 9:221.
- 2.Rodriguez, D. et al. (2010) Biochim. Biophys. Acta 1803:39.
- 3. Bauvois, B. (2012) Biochim, Biophys, Acta 1825;29.
- 4. Heissig, B. et al. (2010) Histol. Histopathol. 25:765.
- 5. Ram, M. et al. (2006) J. Clin. Immunol. 26:299.
- 6. Nagase, H. et al. (2006) Cardiovasc. Res. 69:562.
- 7. Okada, A. et al. (1995) Gene 164:317.
- 8. Xia, Y. et al. (1996) FEBS Lett. 382:285.
- 9. Van Wart, H.E. and H. Birkedal-Hansen (1990) Proc. Natl. Acad. Sci.
- USA 87:5578. 10 Okada, Y. et al. (1992) J. Biol. Chem. 267:21712.
- 11. Collier, I.E. et al. (1992) J. Biol. Chem. 267:6776.
- 12. O'Connell, J.P. et al. (1994) J. Biol. Chem. 269:14967.
- 13. Backstrom, J.R. et al. (1996) J. Neurosci. 16:7910.
- 14. Backstrom, J.R. and Z.A. Tokes (1995) J. Neurochem. 64:1312.

Human MMP-9 Immunoassay

Catalog Number: SEKH-0257

For the quantitative determination of human MMP-9 concentrations in cell culture supernates, serum, and plasma.

For research use only. Not for use in diagnostic procedures.

Country | Company: China | Beijing Solarbio Science & Technology Co.,Ltd Address:NO.85A, Liandong U Valley, Tongzhou District, Beijing, P.R.China.

Tel: 86-10-56371241 Fax: 86-10-56371282 E-mail: service@solarbio.com

TABLE OF CONTENTS

SECTION	PAGE
BACKGROUND	01
PRINCIPLE OF THE ASSAY	01
TECHNICAL HINTS AND LIMITATIONS	02
PRECAUTIONS	02
KIT COMPONENTS& STORAGE CONDITIONS	03
OTHER SUPPLIES REQUIRED BUT NOT SUPPLIED	04
SPECIMEN COLLECTION & STORAGE	04
REAGENTS PREPARATION	04
ASSAY PROCEDURE	06
CALCULATION OF RESULTS	06
PERFORMANCE CHARACTERISTICS	08
REFERENCES	10

LINEARITY: To assess the linearity of the assay, three samples were spiked with high concentrations of MMP-9 in various matrices and diluted with the appropriate Sample Diluent to produce samples with values within the dynamic range of the assay.

Dilution ratio	Recovery(%)	Citrate plasma	Cell culture supernatants
1.2	Average% of Expected	97	105
1.2	Range(%)	88-110	98-112
4.4	Average% of Expected	102	104
1:4	Range(%)	89-113	89-119

Performance Characteristics

SENSITIVITY: The minimum detectable dose was 15 pg/mL.

SPECIFICITY: This assay recognizes both natural and recombinant human MMP-9. The factors listed below were prepared at 100ng/ml in Standard /sample Diluent and assayed for cross-reactivity and no significant cross-reactivity or interference was observed.

Factors assayed for cross-reactivity

Recombinant human	Recombinant mouse	Recombinant porcine
ADAM8	ADAM9	
ADAM10	ADAM10	
ADAM15	ADAM15	
Lipocalin-2/NGAL	Lipocalin-2/NGAL	
MMP-1	MMP-2	
MMP-2	MMP-3	
MMP-3	MMP-9	
MMP-7	TIMP-1	

REPEATABILITY: The coefficient of variation of both intra-assay and inter-assay were less than 10%.

RECOVERY: The recovery of MMP-9 spiked to three different levels in four samples throughout the range of the assay in various matrices was evaluated.

Recovery of MMP-9 in two matrices

Sample Type	Average % of Expected Range(%)	Range(%)
Citrate plasma	97	88-108
Cell culture supernatants	105	92-118


BACKGROUND

MMP-9 is secreted as a glycosylated proenzyme. Activation of the proenzyme involves proteolytic removal of the N-terminal pro region, resulting in the 82 kDa active enzyme. Active human MMP-9 shares 72% and 74% amino acid sequence identity with mouse and rat MMP-9, respectively. In addition to the zinc-binding site, the catalytic domain also contains three contiguous fibronectin type II homology units responsible for binding gelatin. A proline-rich hinge region links the catalytic domain to the C-terminal hemopexin-like domain. In vitro treatment of the proenzyme with 4-aminophenylmercuric acetate (APMA) produces not only the active enzyme but also a C-terminal truncated form with activity comparable to that of the active form). MMP-9 degrades components of the ECM with high specific activity for denatured collagens (gelatin). MMP-9 can also trigger signaling through various transmembrane proteins or inhibit signaling by inducing their shedding from the cell surface (e.g. CD44, E-Cadherin, Integrins, ICAM-1, and IL-2 Rg).

PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. A monoclonal antibody specific for MMP-9 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and anyMMP-9 present is captured by the coated antibody after incubation. Following extensive washing, a biotin-conjugate antibody specific for MMP-9is added to detect the captured MMP-9 protein in sample. For signal development, horseradish peroxidase (HRP)-conjugatedStreptavidinis added, followed by tetramethyl-benzidine (TMB) reagent. Following a wash to remove any unbound combination, and enzyme conjugate is added to the wells. Solution containing sulfuric acid is used to stop color development and the color intensity which is proportional to the quantity of bound protein is measurable at 450nm.

08 www.solarbio.com 400-968-6088 www.solarbio.com 400-968-6088 01



TECHNICAL HINTS AND LIMITATIONS

- 1.This Solarbio ELISA should not be used beyond the expiration data on the kit label.
- 2.To avoid cross-contamination, use a fresh reagent reservoir and pipette tips for each step.
- 3.To ensure accurate results, some details, such as technique, plasticware and water sources should be emphasized.
- 4.A thorough and consistent wash technique is essential for proper assay performance.
- 5.A standard curve should be generated for each set of samples assayed.
- 6.It is recommended that all standards and samples be assayed in duplicate.
- 7.Avoid microbial contamination of reagents and buffers. Buffers containing protein should be made under aseptic conditions and be prepared fresh daily.
- 8. In order to ensure the accuracy of the results, the standard curve should be made every time.

PRECAUTIONS


The Stop Solution suggested for use with this kit is an acid solution. Wear protective gloves, clothing, eye, and face protection. Wash hands thoroughly after handling.

regression analysis. This procedure will produce an adequate but less precise fit of the data. If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.

5. This standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

Typical data using the MMP-9 ELISA

Standared(pg/ml)	OD.	OD.	Average	Corrected
0	0.096	0.091	0.094	
31.25	0.248	0.240	0.244	0.151
62.5	0.310	0.301	0.306	0.212
125	0.479	0.473	0.476	0.383
250	0.844	0.840	0.842	0.749
500	1.400	1.391	1.396	1.302
1000	2.303	2.300	2.302	2.208
2000	3.282	3.109	3.196	3.102



Representative standard curve for MMP-9 ELISA.

DESCRIPTION DESCRIPTION DESCRIPTION

ASSAY PROCEDURE

Prepare all reagents and standards as directed.Wash the plate 3 times before assay.




Add 100µl standard or samples to each well, incubate 90 minutes, 37 °C.



Add 100µl working solution of Biotin-Conjugate anti-human MMP-9 antibody to each well, incubate 60 minutes.37 °C.



Add 100μl working solution of Streptavidin-HRP to each well, incubate 30 minutes,37 °C.



Add 100 μ l Substrate solution to each well, incubate 15 minutes,37 $^{\circ}$ C. Protect from light.



Add 50µl Stop solution to each well. Read at 450nm within 30 minutes.

CALCULATION OF RESULTS

- 1. The standard curve is used to determine the amount of specimens.
- First, average the duplicate readings for each standard, control, and sample. All O.D. values are subtracted by the mean value of blank control before result interpretation.
- 3. Construct a standard curve by reducing the data using computer software capable of generating a four parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph.
- 4.The data may be linearized by plotting the log of the MMP-9 concentrations versus the log of the O.D. and the best fit line can be determined by

KIT COMPONENTS& STORAGE CONDITIONS

PART	SIZE	STORAGE OF OPENED/ RECONSTITUTED MATERIAL
Microwell Plate - antibody coated 96-well Microplate (8 wells ×12 strips)	1 plate	Return unused wells to the foil pouch containing the desiccant pack. Reseal along entire edge of the zip-seal. May be stored for up to 1 month at $2-8$ \mathbb{C}^{**}
Standard - lyophilized,4000 pg/ml upon reconstitution	2 vials	Store at -20°C** for six months
Concentrated Biotin-Conjugated antibody(100X) - 120 ul/vial	1 vial	Store at 2-8°C **for six months
Concentrated Streptavi- din-HRP solution(100X) - 120 ul/vial	1 vial	Store at 2-8°C **for six months
Standard /sample Diluent - 16 ml/vial	1 bottle	Store at 2-8°C **for six months
Biotin-Conjugate antibody Diluent - 16 ml/vial	1 bottle	Store at 2-8°C **for six months
Streptavidin-HRP Diluent - 16 ml/vial	1 bottle	Store at 2-8°C **for six months
Wash Buffer Concentrate (20x) - 30 ml/vial	1 bottle	Store at 2-8°C **for six months
Substrate Solution - 12 ml/vial	1 bottle	Store at 2-8°C **for six months
Stop Solution - 12 ml/vial	1 bottle	Store at 2-8°C **for six months
Plate Cover Seals	4 pieces	

^{**}Provided this is within the expiration date of the kit.

06 www.solarbio.com 400-968-6088 www.solarbio.com 400-968-6088 03

OTHER SUPPLIES REQUIRED BUT NOT SUPPLIED

- 1. Microplate reader capable of measuring absorbance at 450 nm.
- 2. Pipettes and pipette tips.
- 3.Deionized or distilled water.
- 4. Squirt bottle, manifold dispenser, or automated microplate washer.
- 5.500 mL graduated cylinder.
- 6. Human MMP-9 controls (optional; available from Solarbio).

SPECIMEN COLLECTION & STORAGE

Cell Culture Supernates - Centrifuge cell culture media at 1000×g to remove debris. Assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles.


Serum - Use a serum separator tube (SST) and allow samples to clot for 2 hours at room temperature or overnight at 2-8 °C. Centrifuge at approximately for 15 minutes at 1000×g. Assay immediately or aliquot and store samples at \leq -20 °C. Avoid repeated freeze-thaw cycles.

Plasma - Collect plasma using EDTA, heparin, or citrate as an anticoagulant. Centrifuge for 15 minutes at 1000×g within 30 minutes of collection. Assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles.

REAGENTS PREPARATION \

- Temperature returning Bring all kit components and specimen to room temperature (20-25 C) before use.
- 2. Wash Buffer Dilute 30mL of Wash Buffer Concentrate with 570mL of deionized or distilled water to prepare 600mL of Wash Buffer. If crystals have formed in the concentrate Wash Buffer, warm to room temperature and mix gently until the crystals have completely dissolved.

- 3. Standard- Reconstitute the Standard with 1mL of Standard/Sample Diluent. This reconstitution produces a stock solution of 4000pg/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions. Pipette500uL ofStandard/Sample Diluent into 2000pg/ml tube and the remaining tubes. Use the stock solution of 4000pg/mL to produce a 2-folddilution series (below). Mix each tube thoroughly and change pipette tips between each transfer. The 2000pg/mL standard serves as the high standard. The Standard/Sample Diluent serves as the zero standard (0 pg/mL).
 - *If you do not run out of re-melting standard, store it at -20° $\rm C$. Diluted standard shall not be reused.
- 4. Working solution of Biotin-Conjugate anti-human MMP-9 antibody: Make a 1:100 dilution of the concentrated Biotin-Conjugate solution with the Biotin-Conjugate antibody Diluent in a clean plastic tube.
 - *The working solution should be used within one day after dilution.
- Working solution of Streptavidin-HRP: Make a 1:100 dilution of the concentrated Streptavidin-HRP solution with the Streptavidin-HRP Diluent in a clean plastic tube.
 - *The working solution should be used within one day after dilution.



Preparation of MMP-9 standard dilutions

05

04 www.solarbio.com 400-968-6088 www.solarbio.com 400-968-6088