

JM109 Recipient Bacteria

Cat: C1310 Size: 1mL

Storage: Store at -70°C, valid for 1 year.

Introduction:

An amber inhibited F 'recombination defect strain. This strain supports the growth of M13 phage vectors and has modified but unrestricted effects on transfected DNA. F 'in this strain carries lac Z Δ M15, which enables alpha-complementation with the amino end of β -galactosidase encoded in λ ZAP and can be used for blue-white spot screening.

Notes:

- 1. If the biochemical reagents produced by our company are not specially marked, they are basically non-sterile packaging. If they are used in cell experiments, please pre-treat them in advance.
- 2. Once it is prepared into a solution, please pack it separately and store it to avoid product failure caused by repeated freezing and thawing.
- 3. The product information is for reference only, if you have any questions, please call 400-968-6088 for consultation.
- 4. This product is for scientific research only. Do not use for medicine, clinical diagnosis or therapy, food or cosmetics. Do not store in ordinary residential areas.
- 5. For your safety and health, please wear a lab coat and wear disposable gloves and a mask.