

Nucleic acid extraction reagent

Item number: P8440

Specification: 250mL

Storage: Store at room temperature, valid for 2 years.

Product parameters:

CAS: 108-95-2

Appearance (character) : Colorless acicular crystalline solid

Molecular formula: C₆H₆O

Molecular weight: 94.11

Specific gravity: 1.071

Melting point: 40.3°C

Instructions for use:

Put into the water bath, heat to 68°C, do not directly into the 68°C environment, need to slowly heat.

Precautions:

1. Product information is for reference only, if you have any questions, please call 400-968-6088 for consultation.
2. This product is for scientific research only. Do not use for medicine, clinical diagnosis or therapy, food or cosmetics. Do not store in ordinary residential areas.
3. For your safety and health, please wear a lab coat, disposable gloves and a mask.

Related products:

T0250 phenol reagent for DNA extraction

W0250 Reagent for extracting phenol from RNA