

dNTP Mix (25mM each) manual

Item number: PC2500

Grade: PCR grade

Specification: 0.5mL/1mL

Store: -20°C, valid for at least 2 years.

Product Description:

Beijing Solarbio company produced dNTPs (2'-deoxynucleotide-5'-triphosphate) is dNTPs by chemical synthesis, ion exchange chromatography purification, its purity is higher than 98%. 25mM dNTP (2'-deoxynucleotide-5'-triphosphate) Mix contains four deoxynucleotides (dATP, dCTP, dGTP, dTTP), dissolved in 0.6mM Tris-HCl (pH 7.5), each nucleotide concentration is 25mM. 25mM dNTP Mix can be diluted to 10mM or 2.5mM for polymerase chain reaction (PCR), sequencing, filling, incision translation, cDNA synthesis, and TdT tail reaction.

Related products:

PC2200 dNTPs Mix(10mM each)

F0020 Sterile deionized water

PC1100 Taq polymerase

PC1300 Pfu polymerase

M1070 D2000 plus DNA Ladder