

Waymouth MB 752/1 Medium, including Phenol Red and L-glutamine

Cat: C3210

Storage: 2~8°C, away from light.

Validity: 12 months

Product Parameters:

Appearance(character): liquid

Source: Non-animal source

Endotoxin: ≤3 EU/mL

Osmotic pressure: 270-340mOsm/kg

pH: 7.0-7.4

Serum addition: 10% fetal bovine serum Culture condition: 37°C, 5-10% CO₂

Cell type: L929 etc

Notes:

- 1. If the biochemical reagents produced by our company are not specially marked, they are basically non-sterile packaging. If they are used in cell experiments, please pre-treat them in advance.
- 2. Once the solution is prepared, please store it separately to avoid product failure caused by repeated freezing and thawing.
- 3. Product information is for reference only, if you have any questions, please call 400-968-6088 for consultation.
- 4. This product is for scientific use only. Do not use for medicine, clinical diagnosis or treatment, food and cosmetics. Do not store in ordinary residential areas.
- 5. For your safety and health, please wear a lab coat, disposable gloves and a mask.